



tGo

Freedom and mobility for your child
More success with less technology

tGo gait trainer

More success with less technology

- The child is the centre of attention, functionally and visually
- Dynamic saddle supports movement
- Long service life and individual adjustment
- Easy to use even in small spaces



Individual adjustments, simple and straightforward




Flexible pads, individually adjustable



Ergonomically adjustable handle bar



Dynamic saddle that allows the body to rotate



Covered wheel caps reduce risk of injury



Small frame for ease of use



Pads can be extended individually and easy



The friction brake can be adjusted easily with the enclosed tool.



The leg separator plate prevents crossing the legs while walking.



Pull rod



Available in two sizes

Dimensions & weights	Size 1	Size 2
Seat height	25 - 60 cm / 10 - 24"	25 - 80 cm / 10 - 31"
Overall dimensions (L x W)	83 x 54 cm / 33 - 21"	92 x 65 cm / 36 - 26"
Saddle length	20 cm / 8"	23,5 cm / 9"
Handle height	43 - 90 cm / 17 - 35"	45 - 110 cm / 18 - 43"
Pelvic support height	40 - 75 cm / 16 - 30"	40 - 95 cm / 16 - 37"
Wheel size front	15 cm / 6"	15 cm / 6"
Wheel size back	20 cm / 8"	20 cm / 8"
Weight	11 kg / 24.3 lbs	12 kg / 26.5 lbs
Max. capacity	40 kg / 88 lbs	45 kg / 99 lbs

Product & accessories	Size 1	Size 2
tGo gait trainer	1887	1888
Leg separation plate	1889	1890
Divided handle, pair	1893	1893
Vertical hand grips	1895	1896
Pad element (1 piece)	1894	1894
Sit - and - Brake - System	1891	1892
Pull rod	1897	1897