

Bingo rehabilitative chair

catalogue number DRVRX

PRODUCT CHARACTERISTICS:

- wooden rehabilitative chair,
- angle and depth adjustable seat,
- adjustable lateral supports,
- angle adjustable backrest (sitting-horizontal position),
- angle and length adjustable footrests,
- angle adjustable foot platform with heel support and safety belt,
- adjustable abduction wedge,
- height and angle adjustable headrest,
- height and angle adjustable tray,
- positioning vest,
- swivel wheels with lock,
- adjustable headrest,
- height adjustable support of the cervical spine,
- shelf under the seat.

USE:

Rehabilitative chair for disabled children, children diagnosed with nervous system disorder, for children who have difficulty sitting.

Backrest height (with headrest)	46 - 62 cm
Seat depth	25 - 41 cm
Footrest height	25 - 31 cm
Armrest height	68 - 76 cm
Seat width	15 - 35 cm
Backrest width	15 - 27 cm
Chair width	58 cm
Chair length	95 - 150 cm
Chair weight	27 kg
Maximum capacity	35 kg



VITEA
CARE



mdh CE