

www.hoggi.de
HOGGI®

Made in Germany
by HOGGI 

SWINGBO-VTi XL



XL System-wheelchair for adolescents



Tilt from -5° to $+45^{\circ}$

SWINGBO-VTi XL

SWINGBO-VTi XL - Technical data

Seat width	36	40	44
Width			
between seat adapters	35,5 cm	39,5 cm	43,5 cm
Seat depth	34-46 cm	38-50 cm	38-50 cm
Seat height	41-51 cm	41-51 cm	41-51 cm
Seat angle	-5° to +45°	-5° to +45°	-5° to +45°
Backrest height			
angle- and height adjustable backrest	45-54 cm	50-59 cm	54-63 cm
detachable push handle / push bar	43-55 cm	45-57 cm	48-60 cm
Backrest angle			
angle- and height adjustable backrest	75°-115° 80°-120° (in 10° steps)	75°-115° 80°-120° (in 10° steps)	75°-115° 80°-120° (in 10° steps)
detachable push handle / push bar	85°-105° (in 5° steps)	85°-105° (in 5° steps)	85°-105° (in 5° steps)
Lower leg length	20-41 cm	20-41 cm	20-41 cm
Footrest angle	-10° to +10°	-10° to +10°	-10° to +10°
Rear wheel diameter	22" 24"	22" 24"	22" 24"
Front castor diameter	100 mm (4") 125 mm (5") 140 mm (5,5")	100 mm (4") 125 mm (5") 140 mm (5,5")	100 mm (4") 125 mm (5") 140 mm (5,5")
Camber	3°, 6° or 9°	3°, 6° or 9°	3°, 6° or 9°
Load capacity	75 kg	75 kg	75 kg



PLEASE SCAN to see the **SWINGBO-VTi** product video in HD-quality!

HOGGI®

Partnering together to mobilise kids

HOGGI GmbH • Eulerstraße 27 • D-56235 Ransbach-Baumbach • Internet: www.hoggi.de

E-mail: info@hoggi.de • Telefon: (+49) 2623 / 92 499-0 • Telefax: (+49) 2623 / 92 499-99