

www.hoggi.de

Made in Germany
SUPRA
by HOGGI 

SUPRA



Active - Wheelchair for kids & adolescents



• A sporty, rigid and active wheelchair

SUPRA



Nothing is impossible

Axis and backrest independent adjustment **New: SUPRA-S** (frame without abduction)

SUPRA - the light active one with a modular frame system

A sporty, rigid and active wheelchair for kids?

To date, rigid frame active wheelchairs have been reserved for adults or fully grown adolescents.

Once delivered, they have limited adaptability with little or no growth capacity.

To meet the needs of kids, adjustability is the top priority – especially concerning height adjustment.

How can we enable kids to use a sporty and active wheelchair?

The answer is: **SUPRA!**

SUPRA offers a new, modular concept in frame design, creating an active chair that is a suitable and economic choice.

SUPRA - sporty design - economic supply for kids

- The open and rigid frame explicitly represents sportiness.
- The adjustment of axis and the backrest on the seat tubes are independent from one another, enabling ample seat depth adjustment (17 cm) and the ideal center of gravity adjustment (active grade).
- The innovative abduction of the front frame allows for a compact construction of the wheelchair, imparting extra rigidity (for optimal ride characteristics) and stability.
- With a weight of only 8,4 kg for a kids active wheelchair (seat width 26 cm, 100 mm skate wheel), **SUPRA** is a vanguard of active, lightweight wheelchairs.
- It is possible to convert the seat widths from 24 cm to 40 cm, in 2 cm increments, by simply replacing a few transverse members.
- The modular concept provides a sound, economical choice for kids.



Footrest standard

One piece, swingaway footrest

Separated, swingaway footrest

Optional mounting for short lower leg length

HOGGI[®] - Partnering together to mobilise kids



Backrest base in three backrest heights

Height adjustable single push handles

Push bar

High push handles

One-hand push bar

Telescopic tubes for back rest height adjustment (optional)

SUPRA - nothing is impossible

SUPRA is adaptable, offering an extreme combination of seat widths and depths. The modular frame allows for:

- Three frame sizes meet all seat depths, lower leg lengths and seat heights.
- The transverse members (footrest angle, roll bar, rear axle and backrest) can be combined with all three frame sizes.

This results in endless possibilities for seat sizes:

- Seat width 40 cm with a seat depth of 17 cm, or seat width 24 cm with a seat depth of 45 cm.

SUPRA - in detail

- Construction:
 - Rigid frame active wheelchair, aluminium construction, powder coated, adjustable seat depth and backrest.
- 3 size variants:
 - Size 1: SW standard 24 - 30 cm in 2 cm - steps, SD 17 - 34 cm
 - Size 2: SW standard 28 - 34 cm in 2 cm - steps, SD 22 - 39 cm
 - Size 3: SW standard 32 - 40 cm in 2 cm - steps, SD 28 - 45 cm
- Conversion to other seat widths:
- Seat height front:
 - 34 cm - 49 cm
- Seat angle:
 - 0° - 15°
- Backrest angle:
 - 75° - 105°
 - Backrest retractable (optional)
- Foot rest:
 - Footrest standard
 - Footrest with swingaway footboard
 - Footrest with swingaway single footboards
- Brake:
 - Cable brake with hand gear integrated into wheel guard
- Rear wheels:
 - „HOGGI light“ with 1“ air- or PU-tyre 20“, 22“ (20 spokes) and 24“ (24 spokes)
- Front castor wheels:
 - Options with different diameter: 100, 125 and 140 mm
- Weight:
 - From 8,4 kg (**SUPRA** SW 26, 22” drive wheels / 100 mm skate wheels)
- Accessories:
 - Push bar, anti-tip, seat cushion, drum brake system, spoke guard

SUPRA - frame colours

White

Orange

Black

SUPRA

Technical data

Dimensions (cm) & weights (kg)	Size 1	Size 2	Size 3
Seat width			
Standard	24 - 30	28 - 34	32 - 40
Total	24 - 40	24 - 40	24 - 40
Seat width SUPRA-S			
Standard		30 - 34	32 - 40
Total		30 - 40	30 - 40
Seat depth adjustment	17 - 34	22 - 39	28 - 45
Back height adjustment			
Short	20 - 25	20 - 25	20 - 25
Mid	25 - 35	25 - 35	25 - 35
Long	30 - 40	30 - 40	30 - 40
Seat height			
Front	34 - 43	37 - 46	39 - 49
Seat angle	0° - 15°	0° - 15°	0° - 15°
Backrest angle	75° - 105°	75° - 105°	75° - 105°
Lower leg length			
Standard	27 - 41	30 - 44	34 - 48
Rear wheel diameter	20"	20"	20"
	22"	22"	22"
	24"	24"	24"
Front caster diameter	4"	4"	4"
	5"	5"	5"
	5,5"	5,5"	5,5"
Load capacity	75	75	75
Weight	8,4	8,9	9,4

HOGGI®

Partnering together to mobilise kids

HOGGI GmbH • Eulerstraße 27 • D-56235 Ransbach-Baumbach • Internet: www.hoggi.de
E-mail: info@hoggi.de • Telefon: (+49) 2623 / 92 499-0 • Telefax: (+49) 2623 / 92 499-99