

PIPER Comfort

Semi adjustable and foldable standard rehab-buggy



- extended backrest recline
- extended backrest height
- large adjustable footrest
- very compact foldable buggy
- one-piece foldable handlebar
- load capacity up to 50 kg

Size

UNI



PIPER Comfort is the new state of the art comfort rehab-buggy This buggy is determined for short stay usage of children with special needs. Fully foldable undercarriage is produced from lightweight aluminum profiles. Buggy is equipped with seat base with adjustable calf rest, extended reclining backrest, solid wheels of 175/225 mm, big detachable and adjustable footrest, 5-point immo belt with shoulder pads. Backrest is equipped with flexible multi-anchoring system for optimal accessory shoulder alignment. Those features make the products very well suitable for intown activities and public transportation, but as well as for light outdoor activities. The variability of the buggy and the ability to attach a wide range of accessories, ensure fulfilling of majority requirements of user's wide spectrum. The correct usage of the product should be consulted with a seating specialist. Products corresponds to the newest EC standards as well as hygienic, technical and safety requirements.

Standard equipment:

- seat depth non-adjustable
- extended backrest height non-adjustable
- backrest recline adjustable
- footrest height adjustable, detachable
- seat width adjustable by FLEXI-pads set (optional eq.)
- 5-point immo.harness with crotch belt and shoulder padding
- 175/225 mm solid wheels, front with swivel lock
- detachable wheels, forks
- one-piece foldable handle
- basket, transportation lock
- frame colours and cover patterns according to the current offer



folded buggy

Optional equipment:



Flexipads set



Headrest adjustable



Abduction Block



Feet fixation



Canopy Classic



Raincover



Footmuff - Summer



Footmuff - Winter

.... and more according to actual list

Growing functions



Positioning



DESIGN OPTIONS



G06



G07



G12



G13

Order codes:

PIPER Comfort:

PIC34PYYY

	PIPER Comfort
Backrest height [mm]	680
Backrest depth [mm]	290
Seat width [mm]	240-350
Seat to footrest [mm]	250-375
Backrest angle	90° - 135°
Seat angle	26°
Total length [mm]	1080
Total width [mm]	590
Total height/with handle [mm]	1000
Brake	foot oper.
Weight	14,5 kg
Maximum load	50,0 kg
Folded buggy length [mm]	1000
Folded buggy width [mm]	310
Folded buggy height [mm]	300